



Panacea® TRANSPORT DOLLY

for Panacea 3000, 3250 and 3500 beds

Owner's Manual

Please keep and refer to this Owner's Manual.

Thank you for purchasing a Panacea® Transport Dolly from Direct Supply Equipment & Furnishings®. Please read this entire guide carefully and keep it for future reference. This guide will provide you with instructions, warnings, warranty information and other important information about your transport dolly. Share this information with your housekeeping, nursing and maintenance staff to help ensure the transport dolly is cared for properly.

Table of Contents, Special Notes & Notice

Table of Contents

Special Notes.....	2
Notice	2
Warnings	3
Operation.....	4 - 6
Limited Warranty	7
Customer Service.....	8

Special Notes

WARNINGS/CAUTION NOTICES in the manual apply to hazards or unsafe practices which could result in personal injury and/or property damage.

Notice

Check all parts for shipping damage before using. DO NOT USE IF DAMAGED. ONLY USE AUTHENTIC DIRECT SUPPLY REPLACEMENT PARTS. Contact customer service at 1-800-634-7328 and the freight carrier for further instruction.

This product has been engineered to provide reliable operation.

This product has been thoroughly tested and inspected prior to shipment.

Warnings

⚠ WARNING: DO NOT OPERATE THIS EQUIPMENT WITHOUT FIRST READING AND UNDERSTANDING THIS MANUAL. IF YOU ARE UNABLE TO UNDERSTAND THE WARNINGS, CAUTIONS OR INSTRUCTIONS, CONTACT A HEALTHCARE PROFESSIONAL OR A QUALIFIED TECHNICIAN BEFORE USING THIS EQUIPMENT. OTHERWISE INJURY OR DAMAGE MAY RESULT.

⚠ WARNING: NEVER PERMIT ANYONE UNDER THE BED AT ANY TIME.


⚠ WARNING: After installation and before use, make sure all attaching hardware is securely tightened.

⚠ WARNING: The bed is provided with an additional safety feature of a protective earth ground. The electrical cord should only be plugged into a GROUNDED 110/120V outlet. DO NOT ALTER THE PROVIDED PLUG. If the plug does not fit the existing outlet, a licensed electrician should install a correct outlet. Never operate the bed if it has a damaged cord or plug.

⚠ WARNING: Unplug the power cord when performing any maintenance on the bed.

NOTICE: The information in this document is subject to change without notice.

Operation




Required Tools

(2) $\frac{9}{16}$ " wrench


NOTE: Before operating Transport Dolly: headboards/ footboards, siderails and assist rails must be removed. Also, the bed should be raised to its full height.

Attaching Transport Dolly




1. Push the carrier arm brackets into the headboard/ footboard mounting tubes. Lock snap buttons (located on the brackets) into the center holes in the mounting tubes.

NOTE: It is important to push the arm in from both sides to ensure even insertion.




2. Repeat process for the other end of the bed.




3. Using the handset, "lower" bed completely to raise bed casters off floor. Manually stabilize the bed as the casters are raised.
4. Unplug the electrical power cord.

5. Rotate bed ¼ or a turn and align holes between the handle bar and carrier arm.




6. Insert spacer between the handle bar and carrier arm.
7. Insert hex bolt through handle bar, spacer and carrier arm. Tighten jam nut on bolt.
8. Repeat steps 6 & 7 for the other end of the bed. The bed is ready to be transported.




Detaching Transport Dolly


1. Remove jam nut from hex bolt and pull bolt out. Save jam nut and hex bolt.
2. Remove and save spacer.
3. Repeat steps 1 & 2 for the other end of the bed.




4. Rotate the bed so the deck side is up.
5. Remove cable tie from the electrical power cord and handset, and plug the power cord into a compatible receptacle.



Operation *(cont.)*




6. Using the handset, “raise” bed completely to lower bed casters to the floor and raise the dolly wheels off the floor.



7. Press in the snap buttons (locked in the headboard/footboard mounting tubes) and pull out the carrier arm.

NOTE: It is important to pull the arm out from both sides to ensure even removal of brackets.



8. Repeat process for the other end of the bed.

NOTE: The Transport Dolly can be reused, so it is important to reinstall the bolt, spacer and jam nut.

Limited Warranty

We offer to you, as the original purchaser, a warranty for the Panacea Transport Dolly. Our warranty applies for the limited warranty period stated below. If any device or device part listed below is defective in material or workmanship during the applicable limited warranty period, we will repair or replace it at our cost. Please note that the decision to repair or replace a device or device part will be at our discretion. Our warranty applies only if the device is properly maintained by the original purchaser for normal, indoor use and does not cover normal wear and tear, modification of the device, or damage caused by abuse, improper use, failure to maintain, use which exceeds the published device limitations or the combination of any device with another product. In addition, our warranty does not cover fading, colorfastness, stains, spills or exposure to chemicals, odors, heat or light. In certain cases, we may provide you repair or adjustment instructions and/or replacement parts, and ask you to perform a repair or adjustment or replace a defective part.

Our warranty gives you specific legal rights, and you may also have other rights, which vary from state to state. Please note that our limited warranty period begins when we ship the device to you. The limited warranty period and our obligations under the warranty end once you transfer the device to someone else or at the end of the applicable limited warranty period identified below, whichever is earlier.

	Warranty Period (Parts Only)	Anticipated Usable Device Life
Panacea® Transport Dolly	1 year	3 years

Anticipated Usable Device Life is based on normal use with proper maintenance, cleaning and storage. You should still inspect, monitor and care for the device as described in this guide, as the device may need to be replaced sooner than anticipated in particular situations.

DIRECT SUPPLY MANUFACTURING, INC. MAKES NO IMPLIED WARRANTIES OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE; THESE AND ALL OTHER IMPLIED WARRANTIES ARE SPECIFICALLY DISCLAIMED. TO THE FULLEST EXTENT ALLOWED BY LAW, DIRECT SUPPLY MANUFACTURING, INC. WILL NOT BE LIABLE FOR ANY INCIDENTAL, SPECIAL, CONSEQUENTIAL OR PUNITIVE DAMAGES, OR LOST PROFITS THAT MAY RESULT FROM THE DEVICE OR YOUR USE OF OR INABILITY TO USE THE DEVICE EVEN IF ADVISED OF THE POSSIBILITY OF SUCH DAMAGES. OUR TOTAL LIABILITY TO YOU, IF ANY, IS LIMITED TO THE PRICE OF THE DEVICE OR SERVICE GIVING RISE TO YOUR CLAIM. Some states do not allow an exclusion or limitation of incidental or consequential damages, or how long an implied warranty lasts, so the above limitations or exclusions may not apply to you. If implied warranties are not excluded, and to the extent allowed by law, we limit any and all implied warranties to the applicable warranty period identified above. Except for rights under any applicable state law, the remedies provided under this warranty are your sole and exclusive remedy for any breach of our warranty and state the entire limit of our responsibilities.

Customer Service

Our promise to you is that you will have a convenient and easy ordering experience, receive a quality Panacea Transport Dolly and enjoy outrageous customer service. If you have any questions about the bed you have purchased or would like to request warranty service, please contact: **Direct Supply Equipment & Furnishings** at 1-800-634-7328, 6767 N. Industrial Road, Milwaukee, WI 53223, SalesSupport@DirectSupply.com.



DirectSupply.com

Aptura / Equipment & Furnishings / Capital Projects / DSSI / TELS

Distributed by Direct Supply Manufacturing, Inc. 6699 N Industrial Road, Milwaukee, WI 53223.

© 2015 Direct Supply Manufacturing, Inc. All Rights Reserved. Direct Supply®, Panacea®, and all associated logos are the registered trademarks of Direct Supply, Inc., used under license.