

DECONTAMINATION OF CRITICAL OR SEMI-CRITICAL DEVICES PRIOR TO TERMINAL STERILIZATION/HIGH LEVEL DISINFECTION: This product is not to be used as a terminal sterilant/high level disinfectant on any surface or instrument that [1] is introduced directly into the human body, either into or in contact with the bloodstream or normally sterile areas of the body, or, [2] contacts intact mucous membranes but which does not ordinarily penetrate the blood barrier or otherwise enter normally sterile areas of the body. This product may be used to pre-clean or decontaminate critical or semi-critical medical devices prior to sterilization or high-level disinfection.

***THIS PRODUCT KILLS HIV ON PRE-CLEANED ENVIRONMENTAL SURFACES/OBJECTS PREVIOUSLY SOILED WITH BLOOD/BODY FLUIDS** in all settings in which there is an expected likelihood of soiling of hard, non-porous surfaces/objects with blood or body fluids and in which the surfaces/objects likely to be soiled with blood or body fluids can be associated with the potential for transmission of human immunodeficiency virus Type 1 HIV-1.

SPECIAL INSTRUCTIONS FOR CLEANING AND DECONTAMINATING SURFACES/OBJECTS SOILED WITH BLOOD AND BODY FLUIDS WHICH MAY CONTAIN HIV-1 (AIDS virus). PERSONAL PROTECTION: When handling objects soiled with blood or body fluids use disposable latex gloves, gowns, mask and eye coverings.

CLEANING PROCEDURE: Blood and other body fluids must be thoroughly cleaned from surfaces and objects before application of the disinfectant.

CONTACT TIME: Do Not Dilute. Pre-cleaned surfaces must remain visibly wet for five minutes using mop, cloth or sponge to mitigate HIV-1 virus. Use a 5 minute contact time to kill all bacteria, viruses* and fungi listed on the label.

DISPOSAL OF INFECTIOUS MATERIAL: Blood and other fluids should be autoclaved and disposed of according to Federal, State and Local regulations for infectious waste disposal.

Storage and Disposal: Store airtight at room temperature. Non-refillable container. Do not reuse or refill this empty container. Offer empty container for recycling. If recycling is not available, discard in trash.

*Viruses

¹on hard, non-porous surfaces, ²when used according to disinfection directions,

⁴when used as a food contact sanitizer, ⁵100% cotton



INTERCON
For ingredient information,
go to www.rtgclean.com/dfc

Sold by:
Intercon Chemical Co.
1100 Central Ind. Dr.
St. Louis, MO 63110
1-800-325-9218



READY TO GO™

TB Disinfectant Cleaner

Cleans – Disinfects – Tuberculocidal – Fungicidal – Virucidal* – Sanitizer

ACTIVE INGREDIENT:

Citric Acid 0.60%

OTHER INGREDIENTS: 99.40%

TOTAL 100.00%

KEEP OUT OF REACH OF CHILDREN

- Clean and disinfect in one step² • Kills 99.99 % of germs⁴ • Broad Spectrum Hospital Disinfectant
- Effective against Mycobacterium tuberculosis • No rinse, food contact surface sanitizer and cleaner • Sanitizes soft surfaces⁵ in 60 seconds • Suitable for use on hard and/or soft surfaces⁵
- No haze • Dries clear on mirrors and glass • Streak free • Botanically derived ingredients • Free from harsh and abrasive chemicals • Contains no bleach, added dyes, ammonia, phosphates, chlorine bleach, known or suspected endocrine disruptors, glycol ethers, chlorinated compounds
- Non-corrosive • No Rinse required even on children's toys • Note: Not recommended for unpainted wood or natural marble or brass

EPA Reg. No. 34810-35-48211

EPA Est. No. 48211-MO-01

Made in USA



www.epa.gov/pesticides/regulating/labels/design-dfo-pilot.html

NET CONTENTS: 1QT. (32 FL. OZ.)

BACTERICIDAL Bactericidal Activity In the presence of 5 % blood serum, is effective against

5 minutes

Carbapenem-resistant Klebsiella pneumoniae, Enterococcus faecalis Vancomycin Resistant (VRE), Escherichia coli O157:H7 (E. coli), Listeria monocytogenes, Pseudomonas aeruginosa, Salmonella enterica (formerly Salmonella Choleraesuis), Staphylococcus aureus, Staphylococcus aureus - MRSA, Staphylococcus epidermidis - MRSE

VIROCIDAL* In the presence of 5 % blood serum, is effective against

3 minutes

Hepatitis B Virus (HBV), Hepatitis C Virus (HCV), Rhinovirus

5 minutes

Adenovirus Type, Feline Calicivirus (surrogate for Norovirus), Herpes Simplex Virus Type 1, Herpes Simplex Virus Type 2, Human Immunodeficiency Virus Type 1 (HIV-1) (AIDS Virus), Influenza A Virus (flu virus)(H3N2), Respiratory Syncytial Virus (RSV), Rotavirus, Vaccinia Virus

FUNGICIDAL In the presence of 5 % blood serum, is effective against

5 minutes

Candida albicans, Trichophyton interdigitale

TUBERCULOCIDAL At 20°C / 68°F on previously cleaned hard non-porous inanimate surfaces.

5 minutes

Mycobacterium tuberculosis (TB)

FOOD CONTACT SANITIZER 60 seconds

Staphylococcus aureus (Staph), Escherichia coli (E. coli)

SOFT SURFACE SANITIZER 60 seconds**

**Klebsiella aerogenes, **Staphylococcus aureus

DIRECTIONS FOR USE: Please read all directions carefully before using this product. It is a violation of Federal Law to use this product in a manner inconsistent with its labeling. Pre-clean areas with visible soil prior to disinfection.

TO CLEAN: Apply to surface and wipe clean.

TO SANITIZE: For food contact sanitization:

Pre-clean. Wet surfaces. Let stand for 60 seconds on hard, non-porous, food contact surfaces. No rinsing or wiping is required, even on food contact surfaces.

For soft surface sanitization: Pre-clean. Spray this product 6-8 inches from soft surface until wet. Let stand for 60 seconds. Allow to air dry.

TO DISINFECT: Apply the product to hard, non-porous surfaces with sprayer, cloth, mop or brush, or by immersing item to be disinfected. Wet all surfaces thoroughly. The surface must remain visibly wet. Leave for 5 minutes. No rinsing or wiping is required.

FUNGICIDAL DIRECTION: This product is effective against Trichophyton interdigitale (athlete's foot fungus and a cause of ringworm) on inanimate surfaces in locker rooms, dressing rooms, exercise facilities, and shower and bath areas. Follow disinfection directions.

To Spot Clean: Blot or scrub area to be treated with Disinfectant Solution to remove soiling. Wipe dry with clean cloth or towel.

PEEL
HERE

201027-RTBD

Inner Page

On washable hard, nonporous surfaces of: Handrailings, Refrigerators (exterior surface), Respirators¹, Baby/infant care surfaces¹, Sealed Granite, Baked Enamel, Bathing units, Shower floors/walls, Bathroom tubs, Kitchen appliances, Sinks, Beds¹, Kitchen surfaces, Cabinets, Stainless Steel, Chairs, Changing tables, Medical equipment surfaces/devices, Surgical equipment surfaces¹, Children's toys¹, Metal, Chrome, Tables, Computers, Mirrors, Telephones, MRI equipment surfaces, Countertops, Toilets (rims, seats and exterior surfaces), Nonporous athletic mats, Dental equipment surfaces/devices, Desks, Vinyl, Doorknobs, Walls, Exam tables, Exercise equipment¹, Fixtures, Floors, Furniture, Garbage cans, Plastic, Glass, Glazed ceramic tile, Glazed porcelain, PPE (Personal Protective Equipment)¹, as well as other areas/facilities concerned with the hazards of cross contamination

****On the following Soft Surfaces/materials ⁵**

Backpacks, Athletic Equipment, Car seats, Toys, Pillows, Window treatments

AREAS OF USE: Athletic Facilities, Bathrooms, Kitchens, Classrooms, Day care centers, Locker Rooms, Dental offices, Medical Offices, Grocery Stores, Nursing homes, Gyms, Schools, Health Care Facilities, Shopping Centers, Spas, Hospitals, Vehicles, Institutional facilities