

**DECONTAMINATION OF CRITICAL OR SEMI-CRITICAL DEVICES PRIOR TO TERMINAL STERILIZATION/HIGH LEVEL DISINFECTION:** This product is not to be used as a terminal sterilant/high level disinfectant on any surface or instrument that [1] is introduced directly into the human body, either into or in contact with the bloodstream or normally sterile areas of the body, or, [2] contacts intact mucous membranes but which does not ordinarily penetrate the blood barrier or otherwise enter normally sterile areas of the body. This product may be used to preclean or decontaminate critical or semi-critical medical devices prior to sterilization or high-level disinfection.

**\*THIS PRODUCT KILLS HIV ON PRE-CLEANED ENVIRONMENTAL SURFACES/OBJECTS PREVIOUSLY SOILED WITH BLOOD/BODY FLUIDS** in all settings in which there is an expected likelihood of soiling of hard, non-porous surfaces/objects with blood or body fluids and in which the surfaces/objects likely to be soiled with blood or body fluids can be associated with the potential for transmission of human immunodeficiency virus Type 1 HIV-1.

**SPECIAL INSTRUCTIONS FOR CLEANING AND DECONTAMINATING SURFACES/OBJECTS SOILED WITH BLOOD AND BODY FLUIDS WHICH MAY CONTAIN HIV-1 [AIDS virus]. PERSONAL PROTECTION:** When handling objects soiled with blood or body fluids use disposable latex gloves, gowns, mask and eye coverings.

**CLEANING PROCEDURE:** Blood and other body fluids must be thoroughly cleaned from surfaces and objects before application of the disinfectant.

**CONTACT TIME: Do Not Dilute.** Precleaned surfaces must remain visibly wet for five minutes using mop, cloth or sponge to mitigate HIV-1 virus. Use a 5 minute contact time to kill all bacteria, viruses\* and fungi listed on the label.

**DISPOSAL OF INFECTIOUS MATERIAL:** Blood and other fluids should be autoclaved and disposed of according to Federal, State and Local regulations for infectious waste disposal.

**Storage and Disposal:** Store airtight at room temperature. Nonrefillable container. Do not reuse or refill this empty container. Offer empty container for recycling. If recycling is not available, discard in trash.

\*Viruses

<sup>1</sup>on hard, non-porous surfaces

<sup>2</sup>when used according to disinfection directions

<sup>4</sup>when used as a food contact sanitizer

<sup>5</sup>100% cotton



**READY TO GO™**

# TB Disinfectant Cleaner

Cleans – Disinfects – Tuberculocidal – Fungicidal - Virucidal\* - Sanitizer

ACTIVE INGREDIENT:	
Citric Acid .....	0.60%
OTHER INGREDIENTS: .....	99.40%
TOTAL .....	100.00%

## KEEP OUT OF REACH OF CHILDREN

• Clean and disinfect in one step<sup>2</sup> • Kills 99.99 % of germs<sup>4</sup> • Broad Spectrum Hospital Disinfectant  
• Effective against Mycobacterium tuberculosis • No rinse, food contact surface sanitizer and cleaner • Sanitizes soft surfaces<sup>5</sup> in 60 seconds • Suitable for use on hard and/or soft surfaces<sup>5</sup>  
• No haze • Dries clear on mirrors and glass • Streak free • Botanically derived ingredients • Free from harsh and abrasive chemicals • Contains no bleach, added dyes, ammonia, phosphates, chlorine bleach, known or suspected endocrine disruptors, glycol ethers, chlorinated compounds  
• Non-corrosive • No Rinse required even on children's toys • Note: Not recommended for unpainted wood or natural marble or brass

EPA Reg. No. 34810-35-48211  
EPA Est. No. 48211-MO-01

Made in USA



[www.epa.gov/pesticides/regulating/labels/design-dfe-pilot.html](http://www.epa.gov/pesticides/regulating/labels/design-dfe-pilot.html)

NET CONTENTS: 1 GAL ( 3.78L)

**BACTERICIDAL** Bactericidal Activity In the presence of 5 % blood serum, is effective against  
**5 minutes**

Carbapenem-resistant Klebsiella pneumoniae, Enterococcus faecalis Vancomycin Resistant (VRE), Escherichia coli O157:H7 (E. coli), Listeria monocytogenes, Pseudomonas aeruginosa, Salmonella enterica (formerly Salmonella Choleraesuis), Staphylococcus aureus, Staphylococcus aureus - MRSA, Staphylococcus epidermidis - MRSE

**VIROCIDAL\*** In the presence of 5 % blood serum, is effective against  
**3 minutes**

Hepatitis B Virus (HBV), Hepatitis C Virus (HCV), Rhinovirus

**5 minutes**

Adenovirus Type, Feline Calicivirus (surrogate for Norovirus), Herpes Simplex Virus Type 1, Herpes Simplex Virus Type 2, Human Immunodeficiency Virus Type 1 (HIV-1) (AIDS Virus), Influenza A Virus (flu virus)(H3N2), Respiratory Syncytial Virus (RSV), Rotavirus, Vaccinia Virus

**FUNGICIDAL** In the presence of 5 % blood serum, is effective against  
**5 minutes**

Candida albicans, Trichophyton interdigitale

**TUBERCULOCIDAL** At 20°C / 68°F on previously cleaned hard non-porous inanimate surfaces.

**5 minutes**

Mycobacterium tuberculosis (TB)

**FOOD CONTACT SANITIZER 60 seconds**

Staphylococcus aureus (Staph), Escherichia coli (E. coli)

**SOFT SURFACE SANITIZER 60 seconds\*\***

\*\*Klebsiella aerogenes, \*\*Staphylococcus aureus

**DIRECTIONS FOR USE:** Please read all directions carefully before using this product. It is a violation of Federal Law to use this product in a manner inconsistent with its labeling. Preclean areas with visible soil prior to disinfection.

**TO CLEAN:** Apply to surface and wipe clean.

**TO SANITIZE: For food contact sanitization:**

Preclean. Wet surfaces. Let stand for 60 seconds on hard, non-porous, food contact surfaces. No rinsing or wiping is required, even on food contact surfaces.

**For soft surface sanitization:** Preclean. Spray this product 6-8 inches from soft surface until wet. Let stand for 60 seconds. Allow to air dry.

**TO DISINFECT:** Apply the product to hard, non-porous surfaces with sprayer, cloth, mop or brush, or by immersing item to be disinfected. Wet all surfaces thoroughly. The surface must remain visibly wet. Leave for 5 minutes. No rinsing or wiping is required.

**FUNGICIDAL DIRECTION:** This product is effective against Trichophyton interdigitale (athlete's foot fungus and a cause of ringworm) on inanimate surfaces in locker rooms, dressing rooms, exercise facilities, and shower and bath areas. Follow disinfection directions.

**To Spot Clean:** Blot or scrub area to be treated with Disinfectant Solution to remove soiling. Wipe dry with clean cloth or towel.

**On washable hard, nonporous surfaces of:** Handrailings, Refrigerators (exterior surface), Respirators<sup>1</sup>, Baby/infant care surfaces<sup>1</sup>, Sealed Granite, Baked Enamel, Bathing units, Shower floors/walls, Bathroom tubs, Kitchen appliances, Sinks, Beds<sup>1</sup>, Kitchen surfaces, Cabinets, Stainless Steel, Chairs, Changing tables, Medical equipment surfaces/devices, Surgical equipment surfaces<sup>1</sup>, Children's toys<sup>1</sup>, Metal, Chrome, Tables, Computers, Mirrors, Telephones, MRI equipment surfaces, Countertops, Toilets (rims, seats and exterior surfaces), Nonporous athletic mats, Dental equipment surfaces/devices, Desks, Vinyl, Doorknobs, Walls, Exam tables, Exercise equipment<sup>1</sup>, Fixtures, Floors, Furniture, Garbage cans, Plastic, Glass, Glazed ceramic tile, Glazed porcelain, PPE (Personal Protective Equipment)<sup>1</sup>, as well as other areas/facilities concerned with the hazards of cross contamination

**\*\*On the following Soft Surfaces/materials <sup>5</sup>**

Backpacks, Athletic Equipment, Car seats, Toys, Pillows, Window treatments

**AREAS OF USE:** Athletic Facilities, Bathrooms, Kitchens, Classrooms, Day care centers, Locker Rooms, Dental offices, Medical Offices, Grocery Stores, Nursing homes, Gyms, Schools, Health Care Facilities, Shopping Centers, Spas, Hospitals, Vehicles, Institutional facilities



For ingredient information, go to [www.rtgclean.com/dfc](http://www.rtgclean.com/dfc)

Sold by:  
Intercon Chemical Co.  
1100 Central Ind. Dr.  
St. Louis, MO 63110  
1-800-325-9218

201027-RTBD