

Direct Supply® Assist Devices



Multi-Bed Pivoting Assist

- Our most popular assist for Panacea® beds is now compatible with a variety of bed brands and models
- Allows you to standardize on one common assist style, simplifying maintenance, caregiving and resident use
- Unique bracket design includes slotted hole pattern that accommodates some of the most common Long Term Care beds, including models from Panacea, Joerns, Graham-Field and more
- Handle design has a small profile that easily pivots out of the way when not in use
- Also fits beds with 1/2" or 5/8" diameter assist device mounting holes in the head deck (hole spacing supported: 7 1/2", 9 1/4", 10 1/4", 10 1/2", 12 1/2")

#6VY08 Single Assist

#6VY09 Set of 2 Assists



Zero Opening Assist

- A versatile assist device compatible with a variety of bed brands and models
- Sleek, simplistic design features no openings in the handle, minimizing risk of entrapment
- Helps residents reposition themselves in bed and allows for ergonomic hand positioning when egressing from bed
- Slim profile easily pivots out of the way when not in use, and features height adjustment to adapt to resident needs
- Unique bracket design includes slotted hole pattern that accommodates some of the most common Long Term Care beds, including models from Panacea, Joerns, Graham-Field and more
- Also fits beds with 1/2" or 5/8" diameter assist device mounting holes in the head deck (hole spacing supported: 7 1/2", 9 1/4", 10 1/4", 10 1/2", 12 1/2")



#B9G93 Single Assist

#B9G92 Set of 2 Assists



Direct Supply Assist Devices

Tested for compatibility with the following beds:

Bed	Multi-Bed Pivoting Assist	Zero Opening Assist
		
Panacea 3000	X	X
Panacea 3250	X	X
Panacea 3500	X	X
Panacea 1750	X	X
Panacea 1500	X	X
Panacea 8000	X	X
GF Zenith 9100	X	X
Drive P503	X	X
Drive P703	X	X
Medline Alterra 1100	X	X
Medline Alterra Maxx	X	X
NOA Elite Riser	X	X
Panacea 1000	*	*
Panacea 1000-FLEX	*	*
Panacea 2000	**	X
Joerns Care 100	*	*
Joerns EasyCare	**	X

* Pivots toward center of bed only

** Headboard must match deck width

NOTE: Also fits beds with ½" or ⅝" diameter assist device mounting holes in the head deck
(hole spacing supported: 7½", 9¼", 10¼", 10½", 12½")



APTURA / EQUIPMENT & FURNISHINGS / DSSI / TELS

1-800-634-7328 / [DirectSupply.com](https://www.DirectSupply.com)