


Direct Supply® Wheelchair Anti-Rollback Device, Multi-Brand Compatible



#J6K87 – Direct Supply® Wheelchair Anti-Rollback Device



Tool Free Adjustments



Chair can still be folded with the anti-rollback device installed

Benefits

- Intuitive and quick to install, with compatibility on most standard wheelchairs
- Brake arms activate automatically when resident stands, eliminating backward movements to improve safety during ingress & egress. Chair operates normally while resident is seated
- Non-intrusive design allows wheelchair to fold for storage without removing device and does not impede a caregiver's stride when walking behind the chair
- Tool-free adjustment of brake bars, counterbalance bars, and spring force to match many wheelchair designs, with quick-release pins won't loosen or come out of adjustment over time!

#J6K87 Direct Supply® Wheelchair Anti-Rollback Device