

# Direct Supply® Wheelchair Anti-Rollback Device, Multi-Brand Compatible



#J6K87 – Direct Supply® Wheelchair Anti-Rollback Device



*Tool Free Adjustments*



*Chair can still be folded with the anti-rollback device installed*

## Benefits

- Intuitive and quick to install, with compatibility on most standard wheelchairs
- Brake arms activate automatically when resident stands, eliminating backward movements to improve safety during ingress & egress. Chair operates normally while resident is seated
- Non-intrusive design allows wheelchair to fold for storage without removing device and does not impede a caregiver's stride when walking behind the chair
- Tool-free adjustment of brake bars, counterbalance bars, and spring force to match many wheelchair designs, with quick-release pins won't loosen or come out of adjustment over time!
- Patent pending

#J6K87 Direct Supply® Wheelchair Anti-Rollback Device



APTURA / EQUIPMENT & FURNISHINGS / DSSI / TELS

1-800-634-7328 / [DirectSupply.com](https://www.DirectSupply.com)