

NYOrtho™ ZERO-G Heel Pillow ULTRA

The **NYOrtho™ ZERO-G Heel Pillow ULTRA** transforms everyday therapy by providing exceptional support and comfort with a new contoured design that conforms to the body's natural position while off-loading the heel into the air. The ULTRA provides a more relaxing experience, promoting patient position, thereby preventing and providing relief from pressure sores.

RECOMMENDED USE: Individuals unable to avoid long periods of uninterrupted pressure over bony prominences who are at increased risk for the development of necrosis and ulceration. This group of patients typically includes elderly individuals, those who are neurologically impaired, and those who are acutely hospitalized.



ANTIMICROBIAL

No odor, no bacteria: thanks to the built-in hygiene function.



HIGH PERFORMANCE FABRIC

Pressure reducing fabric promotes healing via active evaporation of moisture and cooling at patient pressure points.



WIPE CLEAN

View care label for cleaning instructions.



NON-SLIP

Slip resistant bottom prevents sliding.

- New ULTRA design provides increased comfort and rest while preserving patient position.
- Cushioned design off-loads heels into air, preventing and providing relief from pressure sores.
- Can be positioned in either direction to maintain slightly lifted or lowered position of the heel while providing abduction.
- Convolutd surface securely and comfortably positions both legs, relaxing the spine and decreasing lower back pain.
- Breathable, wipe-clean, low-shear, stretch cover with Antimicrobial protection reduces friction.
- Slip resistant bottom prevents sliding.
- Available strap option secures cushion to bed.
- Meets CA 117-2013 standards.



9902



9901

9902



Model	Description	Dimensions
9901	Single Leg	12 x 17 x 4¾ x 3½ in.
9902	Double Leg	24 x 17 x 4¾ x 3½ in.

Indications for Use

- **Decubitus Ulcer Prevention:** Off-loads the heels to protect against pressure sores and skin breakdown, crucial for at-risk individuals.
- **Aids in Healing Decubitus Ulcers:** Elevates the heels, facilitating easier treatment and re-bandaging of existing wounds, promoting a smoother healing process.
- **Comprehensive Leg Support:** Supports the entire leg including the tibia, ensuring an even distribution of pressure and minimizing the risk of strain.
- **Improves Circulation:** Enhances blood flow back to the heart, beneficial for those with poor circulation or prone to **venous** diseases.
- **Reduces Swelling:** Aids in reducing edema in the legs and feet by facilitating fluid drainage back toward the body's core.
- **Alleviates Varicose Veins Symptoms:** Provides symptomatic relief for individuals suffering from varicose veins by aiding venous return.
- **Enhances Post-Exercise Recovery:** Facilitates the removal of waste products from muscles, aiding in recovery and reducing muscle soreness.
- **Improves Lymphatic Drainage:** Supports the lymphatic system in removing waste and toxins, contributing to a healthier immune system and reducing swelling.
- **Helps with Back Pain:** By allowing for natural spinal alignment, it can relieve pressure on the lower back, offering relief from back pain.

Cleaning, Maintenance and Precautions:

- Clean with a mild detergent and warm water. Do not use abrasive cleaners or solvents, as they may damage the material.
- Allow to air dry completely before using or storing them.
- Store in a cool, dry place away from direct sunlight and heat sources.
- Regularly inspect for signs of wear or damage, and replace them as needed.

Instructions for Use:

- Gently lift the patient's lower legs and place the pillow under them, ensuring that the heels are suspended in the air to off-load pressure without direct contact with the bed. In this position it is easy to treat and bandage the affected are
- Depending on the patient, you may change the desired height/position of the elevated heels by turning the pillow. One side is lower than the other, providing versatility in elevating the heels and optimizing pressure distribution for enhanced comfort and therapeutic benefit.
- When using with optional strap, extend the strap around the bed frame and secure the quick-release buckles.