

# DIRECT SUPPLY® SLING COMPATIBILITY GUIDE



## Direct Supply Sling Compatibility Guide - Floor Lifts

**Direct Supply®** 

Invacare®

### Direct Supply Slings, Multi-Brand Compatible Invacare Hydraulic Reliant 450 Atlas Floor Lift Floor Lift Patient Lift E7793 RPL450-2 9805 Product # **GV270 Description** 9805P RPL450-1 **GV269** #GV258 Full Body Sling, Standard, Extra Small #GV257 Full Body Sling, Standard, Small #GV263 Full-Body Sling, Standard, Medium #GV262 Full-Body Sling, Standard, Large #GV261 Full-Body Sling, Standard, Extra Large #GV255 Full Body Sling, Mesh, Extra Small #GV256 Full Body Sling, Mesh, Small #GV260 Full-Body Sling, Mesh, Medium

Buying compatible slings means no need to spend time searching for the right sling for each

#GV259

#GV254

#GV253

#GV252

#GV251

#GV250

#GV249

#GV248

#GV275

#GV274

#GV273

Full-Body Sling, Mesh, Large

Mesh, Medium

Mesh, Large

Mesh, Extra Large

Divided-Leg Sling, Small

Divided-Leg Sling, Large

Toileting Sling, Small

Toileting Sling, Large

Toileting Sling, Extra Large

Divided-Leg Sling, Medium

Divided Leg Sling, Extra Large

Full Body Sling with Commode Opening,

Full Body Sling with Commode Opening,

Full Body Sling with Commode Opening,



Slings have been tested for compatibility by a third-party accredited test laboratory.

lift in your community; manage your slings inventory with multi-brand compatible slings.

# Direct Supply Sling Compatibility Guide – Floor Lifts cont.

		Joerns®		Graham Field™
	t Supply Slings, -Brand Compatible			
		Hoyer® HPL500 Floor Lift	Hoyer® Presence Floor Lift	Lumex® LF1090 Electric Bariatric Patient Lift
Product #	Description	HPL500-EL-S2 HPL500-S2	HOY-PRESENCE-S	LF1090
#GV258	Full Body Sling, Standard, Extra Small	•	•	•
#GV257	Full Body Sling, Standard, Small	•	•	•
#GV263	Full-Body Sling, Standard, Medium	•	•	•
#GV262	Full-Body Sling, Standard, Large	•	•	•
#GV261	Full-Body Sling, Standard, Extra Large	•	•	•
#GV255	Full Body Sling, Mesh, Extra Small	•	•	•
#GV256	Full Body Sling, Mesh, Small	•	•	•
#GV260	Full-Body Sling, Mesh, Medium	•	•	•
#GV259	Full-Body Sling, Mesh, Large	•	•	•
#GV254	Full Body Sling with Commode Opening, Mesh, Medium	•	•	•
#GV253	Full Body Sling with Commode Opening, Mesh, Large	•	•	•
#GV252	Full Body Sling with Commode Opening, Mesh, Extra Large	•	•	•
#GV251	Divided-Leg Sling, Small	•	•	•
#GV250	Divided-Leg Sling, Medium	•	•	•
#GV249	Divided-Leg Sling, Large	•	•	•
#GV248	Divided Leg Sling, Extra Large	•	•	•
#GV275	Toileting Sling, Small	•	•	•
#GV274	Toileting Sling, Large	•	•	•
#GV273	Toileting Sling, Extra Large	•	•	•

Buying compatible slings means no need to spend time searching for the right sling for each



Slings have been tested for compatibility by a third-party accredited test laboratory.

lift in your community; manage your slings inventory with multi-brand compatible slings.

# Direct Supply Sling Compatibility Guide - Sit-To-Stand Lifts



Buying compatible slings means no need to spend time searching for the right sling for each

### Don't see your lift in this guide?

**Contact your Account Manager for assistance** 



lift in your community; manage your slings inventory with multi-brand compatible slings.

The following spreader bars were used to complete compatibility testing:

Invacare Reliant 450 (RPL450-1, RPL450-2): **1143629** Invacare Reliant 600 (RPL600-1, RPL600-2): **1132034** 

Invacare Hydraulic (9805, 9805P): 1140717

Drive Patient Lift (13240): 13058

Hoyer Hydraulic (HML400) and Hoyer HPL402: **222-600** Hoyer HPL500 (HPL500-S2, HPL500-EL-S2): **31009196** 

Hoyer Presence: 0Y0048

GF Lumex LF1090: **DPL450-SBAR** GF Lumex LF1050: **DPL450-SBAR** 

GF Lumex Hydraulic (LF1030, LF1031, LF1030FP, LF1040, LF1041): LF10406

Medline: MDS400EL: **MDS450CRD**Medline: MDS88200D: **MDS88200DCR** 



1-800-634-7328 DirectSupply.com

**APTURA EQUIPMENT & FURNISHINGS** DSSI TELS

Direct Supply multi-brand slings have been tested with a variety of lifts, as detailed in this guide. When purchasing this sling for use with a non-Direct Supply branded lift, customers must be aware of the following:

1. The use of this sling with a non-Direct Supply branded lift may void the warranty of the lift. Customers are advised to check the warranty period on their lifts before making a sling purchase.

2. State surveyors may require documentation showing that the Direct Supply branded sling is appropriate for use with a non-Direct Supply branded lift. We recommend that you keep a copy of the sling compatibility guide with your

other product manuals and materials.