

Direct Supply Panacea® Air Advance Bariatric Mattresses

High-quality alternating pressure with low air loss for bariatric residents up to 1,000 lbs.

Benefits

- Alternating pressure mode provides dynamic pressure relief
- Low air loss helps wick away heat and moisture to manage maceration
- Each air cell is easily removable for maintenance and replacement
- Alternating pressure cycle times are adjustable from 10, 15, 20 or 25 minutes
- Adjustable firmness levels for resident comfort and clinical efficacy
- Standard firm side perimeters are statically inflated to protect against entrapment and minimize resident falls; 3" raised side perimeters available
- 8"H air cells for maximum resident envelopment and pressure redistribution

Specifications

- May be appropriate through Stage IV wounds with moisture
- 42", 48" or 54"W x 10"H x 80"L
- 750 lb. capacity for 42" and 48"W; 1,000 lb. capacity for 54"W
- 1-year warranty
- UL listed, CA #117

May be appropriate through Stage IV wounds with moisture. Mattress wound stage ratings are general usage guidelines. Resident-specific assessment could alter your particular usage of these mattresses.

SAVE THOUSANDS by owning instead of renting! See details at <u>directsupply.net/mattress</u>

1-800-634-7328

directsupply.com







Direct Supply Panacea® Air Advance Bariatric Mattresses

High-quality alternating pressure with low air loss for bariatric residents up to 1,000 lbs.

Cover

Low-shear, antimicrobial, fluid-resistant cover

- 4-way stretch top cover wicks away heat and moisture, allows for maximum resident immersion
- Bottom cover: heavy-duty vinyl, fluidproof

Concealed zipper to prevent skin tears and protect foam from liquids

- Extra-large flap covers zipper to prevent contact with skin and prevent fluid from seeping in
- Zipper allows for easy cover removal and replacement if necessary

Mattress Body and Side Bolsters

- 20 8"H air cells provide superior envelopment and pressure redistribution
- Each air cell is individually easily removable for maintenance with quick-disconnects
- Air cells are Velcro[®] attached to the top of side perimeters to minimize gaps and resident falls
- Anchor sheet holds air cells securely in the center to prevent them from separating or migrating from their intended location
- Side perimeters and additional bottom layer of 2"H cells are statically filled and remain inflated in the event of a power outage, preventing bottoming out and falls
- CPR quick release valve allows for emergency deflation for superior resident safety

DIFICIT SUPPLY. Equipment & Furnishings





Control Unit

- Provides alternating pressure on a 10, 15, 20 or 25-minute cycle
- 16 lpm low air loss output
- Low pressure light and audible alarm indicates loss of pressure, increasing resident safety
- Allows for use in static low air loss or alternating pressure with low air loss modes
- Adjustable firmness levels for enhanced resident comfort and clinical efficacy
- Max inflate feature automatically adjusts mattress to maximum firmness
- 42"W mattress inflates in 44 minutes

Taking Cost Out of Healthcare®



directsupply.com