



Measuring Knee Length

1. Measure calf circumference at greatest portion to determine size.
2. Measure the distance from bend of knee to bottom of heel (standing, if possible).



Measuring Thigh Length

1. Measure upper-thigh circumference at the buttock fold (if thigh circumference measures greater than 36", select knee length style).
2. Measure calf circumference at greatest portion.
3. Measure distance from base of heel (standing, if possible).



	Thigh Circumference	Calf Circumference	Length	Color		
Thigh Stocking, Short	<25"	Small <12"	<29"	White	#F6996	
Thigh Stocking, Regular		Medium 12" - 15"	29" - 33"		#F6997	
Thigh Stocking, Long		Large 15" - 17½"	>33"		#F6998	
Knee Stocking, Regular	NA	Small <12" Medium 12" - 15" Large 15" - 17 ½" XL 17½" - 20"	Small <16" Medium <17" Large - XL <18"	Beige	#F6999	
				White	#F7001	
				Black	#F7003	
				Nylon	#F7005	
Knee Stocking, Long		Beige	Small >16" Medium >17" Large - XL >18"		#F7000	
					White	#F7002
					Black	#F7004
					Nylon	#F7006
Knee Stocking, Nylon, Regular, XXL - XXXL	XXL 20" - 23" XXXL 23" - 26"		<18"	Nylon	#F7007	
>18"			#F7008			