

DIRECT SUPPLY® SLING COMPATIBILITY GUIDE



Direct Supply Sling Compatibility Guide – Floor Lifts

		Direct Supply®	Invacare®	Drive®		Joerns®
Direct Supply Slings, Multi-Brand Compatible				The state of the s		00000
		Atlas Floor Lift	Reliant 450 Floor Lift	Drive Patient Lift	Drive Levantar Floor Lift	Hoyer® HML400 Hydraulic Lift
Product #	Description	E7793 GV270 GV269	RPL450-2 RPL450-1	13240	FLP500	HML400
#GV258	Full Body Sling, Standard, Extra Small	•	•	•	•	•
#GV257	Full Body Sling, Standard, Small	•	•	•	•	•
#GV263	Full-Body Sling, Standard, Medium	•	•	•	•	•
#GV262	Full-Body Sling, Standard, Large	•	•	•	•	•
#GV261	Full-Body Sling, Standard, Extra Large	•	•	•	•	•
#GV255	Full Body Sling, Mesh, Extra Small	•	•	•	•	•
#GV256	Full Body Sling, Mesh, Small	•	•	•	•	•
#GV260	Full-Body Sling, Mesh, Medium	•	•	•	•	•
#GV259	Full-Body Sling, Mesh, Large	•	•	•	•	•
#GV254	Full Body Sling with Commode Opening, Mesh, Medium	•	•	•	•	•
#GV253	Full Body Sling with Commode Opening, Mesh, Large	•	•	•	•	•
#GV252	Full Body Sling with Commode Opening, Mesh, Extra Large	•	•	•	•	•
#GV251	Divided-Leg Sling, Small	•	•	•	•	•
#GV250	Divided-Leg Sling, Medium	•	•	•	•	•
#GV249	Divided-Leg Sling, Large	•	•	•	•	•
#GV248	Divided Leg Sling, Extra Large	•	•	•	•	•
#GV275	Toileting Sling, Small	•	•	•	•	•
#GV274	Toileting Sling, Large	•	•	•	•	•
#GV273	Toileting Sling, Extra Large	•	•	•	•	•


Buying compatible slings means no need to spend time searching for the right sling for each

Joerns®		Graham Field™			Medline®		Tollos®
					4		
Hoyer® Classics Deluxe Power Floor Lift	Hoyer® HPL500 Floor Lift	Lumex® LF1090 Electric Bariatric Patient Lift	Lumex® Battery Powered Patient Floor Lift	Lumex® Easy Lift Battery Patient Lift	Electric 400-lb. Capacity Patient Lift	Manual Hydraulic Patient Lift	Tollos Ultralift 3000/3510 Floor Lifts
HPL402	HPL500-EL-S2 HPL500-S2	LF1090	LF500	LF1050	MDS400EL	MDS88200D	TUL-600
•	•	•	•	•	•	•	•
•	•	•	•	•	•	•	•
•	•	•	•	•	•	•	•
•	•	•	•	•	•	•	•
•	•	•	•	•	•	•	•
•	•	•	•	•	•	•	•
•	•	•	•	•	•	•	•
•	•	•	•	•	•	•	•
•	•	•	•	•	•	•	•
•	•	•	•	•	•	•	•
•	•	•	•	•	•	•	•
•	•	•	•	•	•	•	•
•	•	•	•	•	•	•	•
•	•	•	•	•	•	•	•
•	•	•	•	•	•	•	•
•	•	•	•	•	•	•	•
•	•	•	•	•	•	•	•
•	•	•	•	•	•	•	•
•	•	•	•	•	•	•	•

Slings have been tested for compatibility by a third-party accredited test laboratory.

lift in your community; manage your slings inventory with multi-brand compatible slings.

Direct Supply Sling Compatibility Guide - Sit-To-Stand Lifts



Buying compatible slings means no need to spend time searching for the right sling for each lift in your community; manage your slings inventory with multi-brand compatible slings.

Don't see your lift in this guide?

Contact your Account Manager for assistance

The following spreader bars were used to complete compatibility testing:

Invacare Reliant 450 (RPL450-1, RPL450-2): 1143629

Drive Patient Lift (13240): 13058

Hoyer Hydraulic (HML400) and Hoyer HPL402: 222-600 Hoyer HPL500 (HPL500-S2, HPL500-EL-S2): 31009196

GF Lumex LF1090: DPL450-SBAR GF Lumex LF1050: DPL450-SBAR Medline: MDS400EL: MDS450CRD Medline: MDS88200D: MDS88200DCR



APTURA

EQUIPMENT & FURNISHINGS

DSSI

TELS

Direct Supply multi-brand slings have been tested with a variety of lifts, as detailed in this guide. When purchasing this sling for use with a non-Direct Supply branded lift, customers must be aware of the following:

The use of this sling with a non-Direct Supply branded lift may void the warranty of the lift. Customers are advised to check the warranty period on their lifts before making a sling purchase.

2. State surveyors may require documentation showing that the Direct Supply branded sling is appropriate for use with a non-Direct Supply branded lift. We recommend that you keep a copy of the sling compatibility guide with your other product