

Westport Recliner

Dimensions:	Power Recliner	Recliner
Overall Height:	39½"	39½"
Overall Depth:	37½"	37½"
Overall Width:	31"	31"
Seat Height:	19½"	19½"
Seat Depth:	21½"	21½"
Seat Width:	20½"	20½"
Arm Height:	24"	24"
Overall Reclined:	66"	66"
Wall Clearance:	17"	17"
Chair Weight:	100 lbs.	85 lbs.
COM Details*:	6 ²² / ₂₅ yds.	6 ²² / ₂₅ yds.
Weight Capacity:	350 lbs.	350 lbs.



Frames:

¾" plywood from managed sources. All joints are interlocking and/or glued, screwed and stapled for added strength and rigidity.

Seat & Back Suspension:

Fully upholstered with a foam covering over a high-tensile-strength steel zig-zag sprung platform.

Seat Cushion:


5" layer of 24 ILD/2.5 density high-resiliency foam wrapped in foam batting.

Back Cushion:


Pillow back with blown foam fiber. Backrest cushions are attached.

Reclining Mechanism:

Power Recliner: Standard power-recline mechanism operated by a corded hand pendant.*



Recliner: Standard 3-position push-back reclining mechanism.



Certifications:

CA117-2013 and CARB II certified. Chair passes BIFMA 5.5X 2014 Seating Standards. CA133 is not available.

Packaging:

Recliner is shipped inside of a box, resting on its feet and wrapped in plastic.


Floor Glides:

Metal.

Limited Warranty:

- 5-year limited warranty on frame and mechanism.
- 1-year limited warranty on foam.
- Fabric manufacturer's warranty on fabric.

Leg Finish Options:



* The estimated fabric requirement is based on 54"W material. The actual amount required may vary depending on pattern, pattern matching, material width or use of multiple fabrics.
* **WARNING:** Do not lean too far forward when the recliner is reclined. This can cause the recliner to tip forward, and could result in injury to the resident and/or caregiver.