

**ENGINEERED
ELECTRONICS, INC.**

**VISUAL PATIENT STATION
900-012R
(Replacement for 900-012)**

FEATURES

- * **AUTOMATIC CALL WHEN CALL CORD IS REMOVED**
- * **DURABLE SINGLE GANG FACEPLATE**
- * **MINIMUM SYSTEM WIRING**



GENERAL DESCRIPTION

The EEI model 900-012R(replacement for 900-012) provides a means of placing a call to the nurses station. This station should be installed in a patient room to accommodate one or two patients which ever the situation requires. All of the solid state call originating circuitry is located on a single gang plate without any subassemblies.

ENGINEERS' SPECIFICATIONS

The EEI model 900-012R are patient stations having all components mounted on a single gang plate.

- | | |
|--|--|
| (1) The equipment described shall be manufactured and assembled entirely in the United States of America | (5) Removal of a call cord either intentionally or by accident shall place a normal call that cannot be reset until the cord is secured in the station. Then the reset switch must be pressed on the station to cancel the call. |
| (2) When a call is placed, a call assurance LED shall illuminate for the patient placing the call. The LED shall continue to glow until reset at the station by momentarily pressing the reset switch. | (6) A normal call must be resetable at the station by momentarily pressing the reset switch located directly beneath the illuminated call assurance LED. |
| (3) The light shall be LED type. | (7) All components must be completely isolated from but mounted to a single gang faceplate requiring two 6/32 machine screws for mounting. |
| (4) Two ¼ inch phone jack receptacles shall be provided to accommodate any type of normally open call cord devices. | |



ASSOCIATED EQUIPMENT


Included

- 2 - 100-023 Dummy Plug
- 2 - 100-195 1-1/4 inch 6/32 countersunk
phillips head machine screw.

Not Included

- 900-800 Oxygen Call Cord 10 ft.
- 900-800-C Oxygen Call Cord 10 ft
- 900-801 Call Cord 6 ft.
- 900-802 Call Cord 10 ft.
- 900-804 Geriatric Presscall Cord 10 ft
- 900-804-C Geriatric Presscall Cord 10 ft
- 900-804-Pad-C Geriatric Pad Call Cord 10 ft
- 900-806 Breath Call Cord 9 ft.
- 900-807 Call Cord 7 ft
- 900-817 Dual Call Cord 10 ft
- 900-819 Call Cord 10 ft W/ Locking Push Button
- 500-188 Pigtail Connector(For Dual outputs)

WIRING DIAGRAM



USEABLE BACK BOXES FOR NEW CONSTRUCTION

Steel City - CDLE, 58371,6060-02 or any single gang outlet box 2 1/2 inches deep with a minimum of 12.5 cubic inch capacity.