

Irvine Recliner

Dimensions:	Power Recliner	Recliner
Overall Height:	40½"	40½"
Overall Depth:	38"	38"
Overall Width:	31"	31"
Seat Height:	20"	20"
Seat Depth:	22"	22"
Seat Width:	19½"	19½"
Arm Height:	24"	24"
Overall Reclined:	67"	67"
Wall Clearance:	18"	18"
Chair Weight:	119 lbs.	104 lbs.
Weight Capacity:	350 lbs.	350 lbs.
COM Details*:	7 yds.	7 yds.



Frames:

¾" plywood from managed sources. All joints are interlocking and/or glued, screwed and stapled for added strength and rigidity.

Seat & Back Suspension:

Fully upholstered with a foam covering over a high-tensile-strength steel zig-zag sprung platform.

Seat Cushion:

6" layer of 24 ILD/1.8 density high-resiliency foam wrapped in foam batting.

Back Cushion:


Pillow back with blown foam fiber. Backrest cushions are attached.

Nailhead Finishes:

Optional.




Frame Finishes:




Reclining Mechanism:

Power Recliner: Standard power-recline mechanism operated by a corded hand pendant.*



Recliner: Standard 3-position push-back reclining mechanism.



Floor Glides:

Metal

Limited Warranty:

5-year limited warranty on frame and mechanism.
1-year limited warranty on foam.
Fabric manufacturer's warranty on fabric.

Packaging:

Chair ships with backrest detached. Backrest is resting face down on top of the seat. Recliner is shipped inside of a box.

Certifications:

CA117-2013 and CARB II certified. Chair passes BIFMA 5.5X 2014 Seating Standards. CA133 is not available.

* The estimated fabric requirement is based on 54"W material. The actual amount required may vary depending on pattern, pattern matching, material width or use of multiple fabrics.
* **WARNING:** Do not lean too far forward when the recliner is reclined. This can cause the recliner to tip forward, and could result in injury to the resident and/or caregiver.