

## Plymouth Recliner

Dimensions:	Power Recliner	Recliner
Overall Height:	42"	42"
Overall Depth:	38"	38"
Overall Width:	29"	29"
Seat Height:	18½"	18½"
Seat Depth:	20"	20"
Seat Width:	21½"	21½"
Arm Height:	26"	26"
Overall Reclined:	67"	67"
Wall Clearance:	6"	6"
Chair Weight:	117 lbs.	95 lbs.
COM Details*:	5½ yds.	5½ yds.
Weight Capacity:	300 lbs.	300 lbs.

### Frames:

¾" plywood from managed sources. All joints are interlocking and/or glued, screwed and stapled for added strength and rigidity.

### Seat & Back Suspension:

Fully upholstered with a foam covering over a high-tensile-strength steel zig-zag sprung platform.

### Seat Cushion:


24-28 ILD/2.5 density high-resiliency foam wrapped in foam batting.

### Back Cushion:


Pillow back with blown foam fiber. Backrest cushions are attached.

### Reclining Mechanism:

**Power Recliner:** Standard power-recline mechanism operated by a corded hand pendant.\*



**Recliner:** Standard 3-position push-back reclining mechanism.



### Certifications:

CA117-2013 and CARB II certified. Chair passes BIFMA 5.5X 2014 Seating Standards. CA133 is not available.

### Packaging:

Chair is shipped inside of a box resting on its feet wrapped in plastic. Backrest is not attached and is resting on the seat.

### Floor Glides:

Metal.

### Limited Warranty:

- 1-year limited warranty on frame and mechanism.
- 1-year limited warranty on foam.
- Fabric manufacturer's warranty on fabric.

\* The estimated fabric requirement is based on 54"W material. The actual amount required may vary depending on pattern, pattern matching, material width or use of multiple fabrics.

\* **WARNING:** Do not lean too far forward when the recliner is reclined. This can cause the recliner to tip forward, and could result in injury to the resident and/or caregiver.