

## Westport Recliner

### Dimensions:

Overall Height:	39½"
Overall Depth:	37½"
Overall Width:	31"
Seat Height:	19½"
Seat Depth:	21½"
Seat Width:	20½"
Arm Height:	24"
Overall Reclined:	66"
Wall Clearance:	17"
Chair Weight:	100 lbs.
COM Details*:	6 <sup>2</sup> / <sub>25</sub> yds.
Weight Capacity:	350 lbs.

### Frames:

7/8" plywood from managed sources. All joints are interlocking and/or glued, screwed and stapled for added strength and rigidity.

### Seat & Back Suspension:

Fully upholstered with a foam covering over a high-tensile-strength steel zig-zag sprung platform.

### Seat Cushion:


5" layer of 24 ILD/2.5 density high-resiliency foam wrapped in foam batting.

### Back Cushion:

Pillow back with blown foam fiber. Backrest cushions are attached.

### Reclining Mechanism:

Standard power-reclining mechanism operated by a corded hand pendant.\*



### Certifications:

CA117-2013 and CARB II certified. Chair passes BIFMA 5.5X 2014 Seating Standards. CA133 is not available.

### Packaging:

Recliner is shipped inside of a box, resting on its feet and wrapped in plastic.


### Floor Glides:

Metal.

### Limited Warranty:

- 5-year limited warranty on frame and mechanism.
- 1-year limited warranty on foam.
- Fabric manufacturer's warranty on fabric.


### Leg Finish Options:



### Available in Westport:



Westport Manual Recliner



Westport Power Recliner

\* The estimated fabric requirement is based on 54"W material. The actual amount required may vary depending on pattern, pattern matching, material width or use of multiple fabrics.  
\* **WARNING:** Do not lean too far forward when the recliner is reclined. This can cause the recliner to tip forward, and could result in injury to the resident and/or caregiver.