



EXTENDED CARE BEDS



We proudly make all of our beds in the USA to the highest safety standards.⁶

- **FDA Recognized Standard:** ANSI/AAMI STD ES60601-1
- **Health Canada Recognized Standard:** CAN/CSA C22.2 No. 60601-1 (IEC 60601-1:2012-Edition 3.1)
- **Includes International Standards:** IEC 60601-1, IEC 60601-1-2, IEC 60601-2-52

All Basic American Extended Care Beds Are

- Certified transport device - passes threshold clearance at ALL heights - ideal for evacuation
- Silver Solutions powder-coating provides antimicrobial surface protection and deters the growth of odor-causing bacteria
- No holes on top surface of tubes that could harbor fluids/bacteria
- Grounded bed with right-angle hospital grade plug for safety
- Adjustable Wallsaver protects the wall with smallest footprint
- Patented Hi-Lo mechanism adds life to the motors
- Sturdy, no-sag, grid sleep surface with cross braces
- Bed provides 360 mobility for cleaning, etc.
- Proudly Made in the USA
- Assist devices available
- Tool-less assembly

Optional
Zenith 7200
Pendant Shown



Adjustable Wallsaver

- Maximum wall protection
- Adjustable clearance lengths
- Standard on all beds



Mattress Retainers

- Adjustable for multiple mattress lengths
- Standard on all beds



Zenith & Matrix Series Bed Standard Features

- Advanced Positioning System (APS) includes One-Button Comfort/Cardiac Chair Position (optional on Matrix 7200).
- Straight-Line Hi-Lo Travel with Patent-Awarded Design Technology
- 'Plug and Play' Digital Electronics enable future feature additions
- Smart Technology
- Industry-leading, energy-saving electronics consume up to 66% less power than transformer electronics when idle
- Sealed ball bearings joints provide continuous, maintenance-free, silent operation
- Intuitive two-pedal locking system - centrally located, easily visible and accessible at all heights and widths
- Tool-less assembly

Zenith & Matrix Series Bed Configurations

- Embedded staff control with safety-inspired, two-stage lockouts enabling caregivers to lock-out features to the resident
- Standard or Advanced Positioning System (APS)
- APS Staff Controls are optional. (APS Staff Control includes One-Button Comfort/Cardiac Chair and Trendelenburg/Reverse Trendelenburg Positions)
- Hand Pendant to control Trendelenburg/Reverse Trendelenburg (*sold separately*)
- Accessories and some options are sold separately



Intuitive Two-Pedal Locking System

- Centrally positioned and easily visible from a distance ensures the caregiver that bed is locked at a glance.
- Bed rolls at any height increasing caregiver safety
- Readily accessible at ALL bed widths and heights (even lowest)
- Select Option at time of purchase

Zenith

Zenith Series Standard Features

- With US Patent awarded time-saver Hi-Lo system
- Tall, Tubular Securely Mounted Head and Footboards create a built-in push bar and reduce breakage

Zenith 9200 SERIES

9200 WITH APS SHOWN WITH
Slide·W·I·D·E® deck option

Zenith 7200 SERIES

7200 & 7200 WITH APS

OPTIONAL WITH THE ZENITH SERIES *Slide·W·I·D·E®* EXPANDABLE/RETRACTABLE SLEEP DECK



Quick-release
ball-detent pins
attached to beds
with lanyards



Wider sleep surfaces reduce falls by residents accustomed to sleeping on common consumer products.²

- Patented design allows for quick adjustment, ideal for moving through doorways to assist in **emergency** evacuations
- Expandable/retractable deck quickly and easily adjusts to 35, 39, or 42 inches with two quick-release ball-detent pins on each side.
- Wide head and footboards available in 42 inch width
- All pins are attached to the bed with lanyards, eliminating lost parts
- Rails and assist devices are also quick-release without tools

² Fragala G, Perry B, Fragala M. Examining bed width as a contributor to risk of falls from bed in long term care. *Annals of Long-Term Care: Clinical Care and Aging*. 2012; 20(6):35-38



Matrix 6200 APS SHOWN
with intuitive pedal-lock option
Raised Head Position

Matrix

Matrix Series Standard Features

- Advanced Positioning System (APS) includes One-Button Comfort/Cardiac Chair Position (optional on Zenith 7200).
- Straight-Line Hi-Lo Travel with Patent-Awarded Design Technology

Matrix 6200 SERIES 6200 WITH APS

Advanced Positioning System (APS)

Advanced positioning allows the bed to go into numerous positions

- Positively influencing clinically driven outcomes through multiple positioning options
- Accommodating post-surgical rehab positioning to facilitate resident comfort
- Reduces caregiver strain and potential injuries versus manually repositioning residents
- Built-In Hi-Lo Lockout Feature (on APS Hand Controls)



Reverse Trendelenburg



Trendelenburg Position



Zenith 9200 APS SHOWN
with Advanced Positioning System, Slide-W-I-D-E®
deck option, and intuitive pedal-lock option
Comfort/Cardiac Chair Position



Liberty SHOWN
with intuitive pedal-lock option
Supine Position



Liberty Bed Standard Features

- Desired features at affordable price point
- Four Locking casters standard

Liberty Bed Available with

- Maximize security with optional one-step red/green pedal lock (*as shown*)
- Optional 4-inch Total-lock casters increase under-bed clearance

Overbed Table Options






Rectangle



Raised Panel



Mission



Suwanee³



Center Crown




Crescent with Cap³



Crescent

Custom Finishes available upon request. Additional cost may incur.

Counter-Rotating Assist Device




Patented Counter-Rotating Assist in Lowered Position



Raised Position

Provide a sturdy and secure hand hold to assist residents in and out of bed yet easily swivels out of the way when desired.

- Meet FDA Entrapment Guidelines
- Patented Counter-Rotating Assist Devices
- Removable without tools
- Move out of the way for lateral transfers
- Pivot away from egress area
- Provide toe clearance in lowest position to prevent staff injury
- Hand pendant mounts facing inward or outward and away from hand grab area



fits Liberty



fits Matrix & Zenith

Pivot Assist Bars



Pivot Assist Device — Rotated Down, Foot



Pivot Assist Device — Up-Right Position




Pivot Assist Device — Rotated Down, Head

Fixed Assist Devices

Provide a sturdy and secure hand hold to assist residents in and out of bed.

- Move out of the way for lateral transfers
- Tool-free assembly and removal



Bed Cart

- Ideal for ease of storage.
- Bed cart allows beds to fit through narrow doorways.
- Designed for One-Person Mounting and maneuverability
- Works with all beds



Side Wall Bumper

- Maximum wall protection
- Works with all Zenith and Matrix Beds

Advanced Positioning Pendant


- Upgrade from Standard to APS on Zenith 7200
- One-Button Comfort/Cardiac Chair Position
- Hi-Lo Lockout Feature



USB Power Supply

- Charges most devices
- Works with all beds





Board Mounted IV Socket

- Mounts with 2 screws
- Easy to install
- IV pole sold separately
- Works with all beds




Plastic Pan Deck Kit

- Fits 6200, 7200, 9200 Series



Edema Foot Ratchet Kit

- Elevated foot rest.
- Standard on Zenith 9200
- Optional on Matrix 6200 and Zenith 7200




Surface Length Extension Kit

- Extends standard length bed frame from 84" to 88". Complete with solid pan extension, mattress guard, and installation hardware.
- Works with all beds
- Easy to install



UnderBed Light

- Reduces fall risks
- Eliminates bright light intrusion
- Multi-level light intensity
- Works with all beds



Pin-on Wide Kit

- Fits 35 inch width 6200, 7200, 9200
- Available in 39" or 42" width kits
- Anytime add-on



Embedded Staff Control

- Standard Functionality and Lockouts, or Advanced Positioning (APS).
- Angled-Up Display for Improved Visibility
- Two-Stage lockouts deter unwanted changes
- Mix or Match with Plug & Play design
- Upgrade anytime



Trapeze and Adapter

- Pivots out of the way when not in use
- Works with all Zenith and Matrix Beds
- Rigid Design



 Accessories

Accessories Sold Separately

Kensington
Maple

Solar Oak

Fonthill Pear

Wild Cherry

Windsor
Mahogany

Shaker Cherry

Talas Cherry

Aria

Cocobala

Bed Specifications:

	Zenith 9200	Zenith 7200	Matrix 6200	Liberty
Adjustable Patient Surface Length:	80-84 inch	80-84 inch	80-84 inch	76-80 inch
Welded Sleep Deck Width:	35 inch	35 inch	35 or 42 inch	35 inch
Height Travel Range:	7-30 inch	7.95-30 inch	8.95-27 inch	9.95-24.4 inch
Safe Work Load: <i>includes resident, support surface, head/footboards, and accessories</i>	600 lb	550 lb	450 lb	450 lb
Maximum Patient Weight: EVENLY DISTRIBUTED	550 lb	500 lb	400 lb	400 lb
Actuators:	4 DC	4 DC	4 DC	3 DC

Options & Features:

Chair Position (APS) :				—
Advanced Positioning Pendant:				—
Pendant Attaches to Either Side of Bed:				—
Trendelenburg/Reverse Trendelenburg Positioning: via APS staff control				—
Slide-W-I-D-E™ Option Available: between 35 - 42 inches			—	—
Patented Super-Efficient Transmission: less strain on the motors		—	—	—
Hi-Low Actuators are Unloaded in Low Position:				—
Staff Control: ¹				—
Upgrade Hand Controls and/or Staff Controls Anytime: interchangeable, mix and/or match standard and/or APS controls				—
Splash Proof Electronics:				—
Intuitive Pedal Lock or 4 Locking Casters:				 ¹
Bed Locks and Rolls at Any Height: reducing core-giver back strain				
Silent Operating, Sealed Ball Bearings: Maintenance-Free				—
Taller Casters: approximately 7 inches of clearance	—	—	—	
Head and Footboard Mounting:	Tall Tubular	Low Profile	Low Profile	Low Profile
Built-in Push Bar and Breakage Reduction:		—	—	—
Maximum Clearance: over bed, tables and lifts				51 inch ⁴
Can be Used as Transport Devices:				—
FDA Guideline Compliance: for Bed Entrapment Risk Reductions				—
Silver Solutions Powder Coating: process provides antimicrobial surface protection				—
Made in the USA by Skilled Craftsmen: built to withstand trouble-free use for many years				—

Accessories:

Edema Foot Ratchet: (elevates the foot end of the bed)				—
Plastic Pan Deck:				—
Backup Battery: (external or mounted onboard)				—
Extension Kit: (converts from 84 to 88 inch length)			—	—
Pin-on Wide Kit:				—

Warranty:

Main Frame & Welds	15 years	15 years	15 years	10 years
Control Box & Actuator Motors	5 years	5 years	4 years	2 years
Hand Control Pendant & Cabling	3 years	3 years	3 years	1 year
Staff Control & Cabling	3 years	3 years	3 years	—
Head & Footboards	1 year	1 year	1 year	1 year

Product Page Link

			
Standard	Available	Optional	Not Applicable/Available

¹ anytime add-on or with purchase³ boards made with 3DL material - all styles available in 3DL - ask for details — custom finishes available upon request (additional cost may incur)⁴ fixed base⁵ "Hospital Bed System Dimensional and Assessment Guidance to Reduce Entrapment" details: www.fda.gov⁶ preserves the cured coating during the service life of the bed & inhibits the growth of organisms that can detract from the appearance and coating functionality. inhibits the growth of bacterial odors and resists stains caused by mildew and/or mold

GF Health Products, Inc. makes no independent representations or empirical statements about Silver Solutions. This information is reprinted from: TCI Powder Coatings, EPA Guidance for EPA Registration 40 CFR 152.25

GF Health Products, Inc. is an
ISO 13485:2016 / MDSAP Certified Company Color variations may occur due to monitor or printer settings.

+1 770.368.4700

Information contained herein is subject to change without notice.
The most current and complete product information can be found on our website.www.grahamfield.com© 2017, GF Health Products, Inc. All Rights Reserved.
Basic American Medical Products is a trademark of GF Health Products, Inc.


GF700011RevH23