Panacea® Air Pro Mattress





The Panacea® Air Pro alternating-pressure and low air loss mattress is designed to assist in the prevention and treatment of pressure ulcers for Long Term Care residents. The resident support surface fits most standard healthcare beds and consists of 20 air cells that run laterally underneath the entire body and are contained within durable air perimeters. The air flotation system aids in pressure redistribution and relief, and helps maintain low interface pressures throughout the surface.

The Panacea Air Pro is a pump-based system that provides alternating-pressure and low air loss therapy. The inflation levels are designed to maintain low interface pressures throughout the mattress, to redistribute peak interface pressures during the alternating cycle and to further relieve peak pressure points as the cycle alternates. The true low air loss feature helps reduce maceration and further skin breakdown by keeping the skin surface cool and dry.

	SIZE		CONFIGURATION				
	Width	Length	Flat	Raised Perimeter	Heel Slope	Raised Perimeter & Heel Slope	
PANACEA® AIR PRO	35"	75"	#A9160	#A9164	#A9162	#A9166	
		80"	#A9161	#A9165	#A9163	#A9167	

See	DirectSupply n	at for	hariatric	99712

Specifications				
Standard Dimensions	35"W; 11"H; 76" or 80"L			
Weight Capacity	500 lbs.			
Therapy Style	Alternating Pressure, Limited Low Air Loss			
Firm Perimeters	Yes			
Raised Perimeters	Yes, Optional			
Heel Slope	Yes, Optional			
CPR Release	Yes			
Safety Base [†]	Yes			
Air Cell Height	8"			
InflationTime	Less than 60 minutes			
CoverType	4-way stretch (Integrated Low Air Loss)			
Limited Warranty	Mattress – 2 Years; Pump – 18 Months			

[†]Prevention against bottoming out during loss of power or air cell puncture

These pressure redistribution support surfaces may be appropriate for use as part of an overall care plan to help prevent and treat decubitus ulcers. Resident-specific assessment could alter your particular usage of these mattresses.



DirectSupply.com