

Panacea® Air Ease Mattress



#2TPDF Panacea® Air Ease Mattress



Benefits

- The pump offers a 10-minute cycle time; adjustable comfort settings allow pressure to be customized by resident's weight.
- The static feature offers constant low-pressure support and an internal pressure adjustment.
- Audio and visual alarm indicates low pressure while cell-on-cell design slows deflation time in the event of a power outage.

#2TPDF Panacea® Air Ease Mattress

Specifications	
Standard Dimensions	35"W x 8"H x 79"L
Weight Capacity	485 lbs.
Therapy Style	Alternating Pressure, Limited Low Air Loss
Firm Perimeters	No
Raised Perimeters	No
Heel Slope	No
CPR Release	Yes
Safety Base [†]	No
Air Cell Height	8" (4" cell-on-cell design)
Inflation Time	Less than 45 minutes
Cover Type	2-way stretch
Limited Warranty	Mattress – 1 Year; Pump – 2 Years; Cover – 1 Year

[†]Prevention against bottoming out during loss of power or air cell puncture.

These pressure redistribution support surfaces may be appropriate for use as part of an overall care plan to help prevent and treat pressure injuries. Resident-specific assessment could alter your particular usage of these mattresses.

