

Panacea® Air Overlay System

Bubble pad design for superior treatment



#36BCK Air Overlay System with End Flaps and Adjustable Pump

#36BCK Air Overlay System – end flaps tuck under mattress for secure fit

Benefits

- Bubble pad design provides more touch points that cycle to reduce pressure
- Adjustable pump allows you to customize the amount of pressure
- End flaps tuck under mattresses to secure pad to mattress
- Small, compact design
- 250 lb. weight capacity
- UL listed, with hospital-grade plugs

#36BCK Panacea® Air Overlay System, with End Flaps and Adjustable Pump

Accessories

#36BCM Replacement Adjustable Pump for Panacea® Air Overlay

#36BCP Replacement Panacea® Air Overlay, Pad Only

#36BCJ Pad Repair Kit for Panacea® Air Overlay

Specifications	
Size	35"W x 79"L x 2.5"H
End Flaps	Yes
Adjustable Pump	Yes
Weight Capacity	250 lbs.
Bubble Layers	Single
Air Flow	6 - 8 lpm
Cycle Time	10 minutes
Alternating Pressure	Yes
Pressure Range	30 - 120 mm Hg
Meets CA #117, UL 1431	Yes
Pump Weight	3.3 lbs.
Pump Dimensions	10.75"W x 4.75"L x 4"H
Overlay Material	Polyvinyl chloride
Foldaway Pump Hanger	Yes



APTURA / EQUIPMENT & FURNISHINGS / DSSI / TELS

1-800-634-7328 | DirectSupply.com